
# PentaView™

LCD DIGITAL MICROSCOPE

#44348

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# INTRODUCTION

Thank you for purchasing the Celestron PentaView™ LCD Digital Microscope with a 4.3-inch touch screen monitor. Your microscope is a precision optical instrument, made of the highest quality materials. It is designed to give you a lifetime of enjoyment with minimal maintenance.

Before attempting to use your microscope, please read through the instructions to familiarize yourself with its functions and operations. See the microscope diagram to locate the parts discussed in this manual.

This microscope provides high powers from 40x up to 600x, and up to 2400x with the digital zoom. This microscope is mainly suited for examining specimen slides of yeasts and molds, cultures, plant and animal parts, fibers, bacteria, and more.

Unlike traditional microscopes, the Celestron PentaView™ LCD Digital Microscope uses a 5MP camera sensor in lieu of an eyepiece. You will view specimens on the LCD screen. In addition, the sensor allows you to take snapshots or short videos. Plus, you can view on most TV screens with the AV/TV cable.

The final sections of this manual provide simple care, maintenance and troubleshooting tips to ensure that your microscope provides you with years of quality performance and enjoyment.

Fig 1

## FEATURES

- Built-in Digital Camera
- 4x Objective Lens
- 10x Objective Lens
- 20x Objective Lens
- 40x Objective Lens
- 60x Objective Lens (Included Separately)
- LCD Touch Screen Monitor
- Illuminator
- 6 Position Filter Wheel/Diaphragm
- SD Card Slot

## INCLUDED ACCESSORIES

- 10 Prepared Slides
- Hard Carrying Case
- Dust Cover
- AC Adapter
- 8 GB SD Card
- USB Cable 2.0
- AV/TV Cable
- Touch Pen

## SPECIFICATIONS

|                          |                                                                                                                         |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------|
| <b>Stage</b>             | Mechanical Stage 3.5 in x 3.5 in (88 mm x 88 mm)                                                                        |
| <b>Digital Camera</b>    | 5 MP 1/2.5" CMOS; 10x Magnification in lieu of an eyepiece                                                              |
| <b>LCD Monitor</b>       | 4.3 in. (109 mm) TFT LCD touch screen with 4x digital zoom – Touch digital TFT display<br>Resolution – 480 x 272 pixels |
| <b>Focuser</b>           | Dual – Coarse/Fine                                                                                                      |
| <b>Objectives</b>        | Achromatic, glass type – 4x, 10x 20x, 40x and 60x                                                                       |
| <b>Memory</b>            | 8 GB SD card (approximately 2200+ snapshots @ 5 MP)                                                                     |
| <b>LCD Rotation</b>      | 180° – 90° left / 90° right                                                                                             |
| <b>Filter Wheel</b>      | Red/Green/Blue/1 mm hole/3 mm hole/6 mm hole (default)                                                                  |
| <b>Nosepiece</b>         | Quad with click stop                                                                                                    |
| <b>Illuminator</b>       | Built-in electric – LED 6 Volt and 6 Watt                                                                               |
| <b>Condenser</b>         | N.A. 0.65                                                                                                               |
| <b>AC Adapter</b>        | Input Universal 100 to 240 Volt 50/60HZ                                                                                 |
| <b>Batteries</b>         | User supplied 4 AA – (optional - up to 3 hours usage)                                                                   |
| <b>Weight/Dimensions</b> | 67 oz/1.9 kg 7.0" (178 mm) x 5.5 in (140 mm) x 13.0 in (330 mm)                                                         |

## MAGNIFICATION (POWER)

Use the following table to determine the magnification of the different objective lenses in combination with your microscope using the normal mode of the digital image on the LCD screen and using the digital zoom feature.


|                                      |      |      |      |       |       |
|--------------------------------------|------|------|------|-------|-------|
| Objective Lens                       | 4x   | 10x  | 20x  | 40x   | 60x   |
| Digital Image                        | 40x  | 100x | 200x | 400x  | 600x  |
| Maximum with 4x digital zoom feature | 160x | 400x | 800x | 1600x | 2400x |

## WORKING DISTANCE

|                              |      |     |     |     |     |
|------------------------------|------|-----|-----|-----|-----|
| Objective Lens               | 4x   | 10x | 20x | 40x | 60x |
| Working Distance (mm) - est. | 35.3 | 7.8 | 1.9 | 0.7 | 0.2 |

## TOUCH PEN

You can keep the touch pen handy in the receptacle behind the Power ON/OFF button.

## SETTING UP YOUR MICROSCOPE

1. Carefully remove the microscope and other parts from the carton and set them on a table, desk, or other flat surface.
2. Remove the plastic cover from the microscope.
3. Plug the small cable from the AC adapter into the socket on the back of the base (see Figure 2).

Fig 2

4. Insert the plug you need (see note below) into the AC adapter and then plug the adapter into the proper power source.

**NOTE:** The AC Adapter supplied with this microscope has an interchangeable plug system that can work in the USA, UK, Europe, Australia, and many other parts of the world.

You can easily change plugs for your particular power requirement:

- A. You may have one of two different types of adapters. One is shown in Figure 3A and another one is somewhat similar shown in Figure 3B. Each has a different means of changing plugs.
- B. For the Figure 3A type, push down on the small button labeled "PUSH" and hold it down while at the same time holding the prongs of the plug. Rotate it slightly counterclockwise and pull up to remove it. Then release the "PUSH" button down.

Take the plug style that you want to use and center it over the AC adapter and rotate it until it falls into the opening. Rotate it clockwise until you hear a click, which lets you know the plug is installed correctly.

Fig 3A

- C. For the Figure 3B type, pull the plug up and away while holding the base of the AC adapter. The plug will come off. Take the plug style that you want to use and align the two holes in the bottom of the plug with the two pins in the AC adapter. Then, push the plug into the AC adapter all the way.

Fig 3B

**Battery Operation** – You can use your microscope without AC power if you choose. This gives you the freedom to operate the microscope outdoors or indoors anywhere you please. You'll need 4 AA batteries (user supplied). Open the battery door on the bottom of the microscope and insert the batteries according to the battery polarity shown in the battery compartment (Figure 4). After the batteries are installed, close the battery door. Battery life will typically be three hours.

Fig 4

Product design and specifications are subject to change without prior notification. This product is designed and intended for use by those 14 years of age and older.

© 2019 Celestron • All rights reserved • celestron.com/pages/technical-support • 2835 Columbia Street • Torrance, CA 90503 USA • 12-19  
celestron.com/support/warranties

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